

## ParaSure DentaPro

- Unique & unrivalled weekly residual protein test
- Cost effective
- Accurate results within 10 seconds
- No other test more sensitive without incubation
- HTM 01 05 compliant
- Easy to use
- Clear colour change when protein detected
- Enhanced accuracy
- Shelf life 2 years



## ParaScope DentaPro

- The ParaSure DentaPro protein test is a rapid test that has the ability to detect residual proteins left behind on the surfaces of hard to clean dental and surgical instruments. An easy to read colour change gives you your final result within 10 seconds. DentaPro is suitable for use with washer disinfectors, ultrasonics and even manual cleaning.

- This protein test is based on a dye-binding solution used in clinical chemistry and can detect protein residues within 1µg sensitivity, therefore incubation is not required.

- DentaPro meets and exceeds current guidelines for protein testing, as published in HTM 0105 (pg 60 14.4)

Order Details

**PS1740**

ParaSure DentaPro

Box of 10

## ParaSure DentaPro

TECHNICAL SPECIFICATIONS	
<b>Name:</b>	DentaPro
<b>Conformity:</b>	BS EN ISO 15883, HTM 2030
<b>Time:</b>	Results within 10 secs
<b>Storage conditions:</b>	Refridgerated (4°C)
<b>Lot no:</b>	yes
<b>Exp date:</b>	yes
<b>Correct symbols:</b>	yes in accordance with the standard
<b>Indicator:</b>	Diagnostic dye-binding vial
<b>Expiry date:</b>	24 months refridgerated, 6 months non-refridgerated
<b>Packaging:</b>	10 tests per pack including swabs + 1 protein control test

Quick and easy to use; to determine if protein is present on a surface swab the instrument surface, place the swab in the colour changing solution and swirl for 5 seconds. The swab can then be discarded. If protein is present the reagent will turn blue within 10 seconds. The colour change provides an easy to read result; the higher the level of protein contamination the darker the colour changes to blue.



Address: Squirrels | Kithurst Park |  
Storrington | West Sussex | RH20 4JH | England |  
Tel: +44 (0) 1903 741390

Warehouse Address: Units B10/B11  
Oakendene Estate | Bolney Road | Cowfold |  
West Sussex | RH13 8AZ | England