



PARASURE

STOP!

be sure...be ParaSure

Type 4 Multi-Parameter Indicators

Unprocessed indicator Large square is pale pink



Successfully processed indicator

All lines shows evidence of black.

PASS

Items in tray/pack can be safely used



Unsuccessfully processed indicator

one or more lines shows evidence of pink.

FAIL

Do NOT use items in tray/pack



Unsuccessfully processed indicator

one or more lines shows evidence of pink.

FAIL

Do NOT use items in tray/pack



Unsuccessfully processed indicator

one or more lines shows evidence of pink.

FAIL

Do NOT use items in tray/pack

PS1110 Type 4 Multi-Parameter Indicators



WWW.PARASURE.CO.UK

Registered Address:

ParaSure Limited Units B10/11 Oakendene Estate
Bolney Rd Cowfold West Sussex RH13 8AZ England



WWW.PARASUREORGANICS.COM

+44 (0)1403 864777



PARASURE

STOP!

be sure...be ParaSure

Type 4 Multi-Parameter Indicators

Unprocessed indicator Large square is pale pink



Successfully processed indicator

Large square is as dark or darker than small square

PASS

Items in tray/pack can be safely used



Unsuccessfully processed indicator

Large square shows evidence of pink

FAIL

Do NOT use items in tray/pack



Unsuccessfully processed indicator

Large square shows evidence of pink

FAIL

Do NOT use items in tray/pack



Unsuccessfully processed indicator

Large square shows evidence of pink

FAIL

Do NOT use items in tray/pack

PS1110 Type 4 Multi-Parameter Indicators



WWW.PARASURE.CO.UK

Registered Address:

ParaSure Limited Units B10/11 Oakendene Estate
Bolney Rd Cowfold West Sussex RH13 8AZ England



WWW.PARASUREORGANICS.COM

+44 (0)1403 864777