

PARASURE

PS1610

Standard Basket PCD Process Challenge Device



*For monitoring the cleaning efficacy of
Washer Disinfectors to EN ISO 15883-5:2005*

*To be used in conjunction with
ParaSure Cleaning Efficacy Indicators*

For further information about any of our products, please visit:

www.parasure.co.uk

ParaSure Limited
Units B10/11 Oakendene Estate | Bolney Road | Cowfold | West Sussex | RH13 8AZ | England
Tel: +44 (0) 1903 741390 Email: info@parasure.co.uk

PARASURE

Standard Basket PCD

PS1610